

# Irrigation Sets

Fresenius Kabi offers a range of sterile irrigation sets to complement our **soluflex™** Irrigation Bag portfolio and to support health care providers in delivering a high standard of care.



Fresenius Kabi irrigation sets are available in both 1-spike and 2-spike configurations, each designed to offer distinct features depending on the specific set. These features may include:

- Bag spike with drip chamber
- Roller clamps
- Urology connector
- Backcheck valve
- Pinch clamps
- Luer lock connection
- Drain hose with Luer lock
- Silicone para tube (allows connection between the urology connector and other components)



Contains DEHP



Not made with DEHP


All sets are not made with natural rubber latex and contain PVC

## Single Lead Urology Sets

### TRANSFER SET 1 D

Drip chamber with spike, roller clamp, urology connector, and optional para tube.


Product Code	GTIN	Length (inch)	Shelf Life (months)	Units / case
2883101	04250273756081	70.9	60	30



### TRANSFER SET 6 D

Drip chamber with spike, roller clamp, backcheck valve, urology connector, and optional para tube.

Product Code	GTIN	Length (inch)	Shelf Life (months)	Units / case
2884201	04250273755909	70.9	60	30




## Single Lead Irrigation Sets

### TRANSFER SET 4 D/X (double sterile)

**First tube:** Drip chamber with spike, roller clamp, and positive Luer lock connector.

**Second tube:** Drainage tube with positive Luer lock connector, and glued para tube.


Product Code	GTIN	Length (inch)	Shelf Life (months)	Units / case
2848021	04250273755749	67, 51.2	60	15



### TRANSFER SET 2 D/X (double sterile)

Drip chamber with spike, roller clamp, and positive Luer lock connector.

Product Code	GTIN	Length (inch)	Shelf Life (months)	Units / case
2877201	04250273755787	130	60	15




All measurements are approximate.

## Dual Lead Urology Sets

### TWIN TRANSFER SET 1 D

2 spikes, 2 roller clamps, backcheck valve, urology connector, and optional para tube.


Product Code	GTIN	Length (inch)	Shelf Life (months)	Units / case
2894201	04250273738629	128	36	15



### TWIN TRANSFER SET 1 D/X (double sterile)

2 spikes, 2 roller clamps, backcheck valve, urology connector, and optional para tube.


Product Code	GTIN	Length (inch)	Shelf Life (months)	Units / case
2877001	04250273756005	126	60	10



### TWIN TRANSFER SET 2 D

2 drip chambers with spikes, 2 roller clamps, urology connector, and optional para tube.


Product Code	GTIN	Length (inch)	Shelf Life (months)	Units / case
2890101	04250273756029	76.8	36	15



### TWIN TRANSFER SET 9 D

2 spikes, 2 roller clamps, inline drip chamber, 1 pinch clamp, urology connector, and optional para tube.

Product Code	GTIN	Length (inch)	Shelf Life (months)	Units / case
2890701	04250273756067	98.4	36	15




## Dual Lead Irrigation Sets

### TWIN TRANSFER SET 2 D/X (double sterile)

**First tube:** 2 drip chambers with spikes, 2 roller clamps, and positive Luer lock connector.

**Second tube:** Drainage tube with positive Luer lock connector, roller clamp, and glued para tube.

Product Code	GTIN	Length (inch)	Shelf Life (months)	Units / case
2848011	04250273755725	114.2, 122	36	8



All measurements are approximate.

# Together, we can care better

Get to know our full line of IV Therapy products



**IVENIX**  
INFUSION SYSTEM

Engineered from the ground up to leverage the power of the digital age, Ivenix system is designed with safety in mind.



**agilia**  
Connect

Agilia Connect integrated system of pumps, software, sets, accessories and services.



**freeflex**<sup>®</sup>  
IV Bags & Solutions

A growing portfolio of IV solutions in **freeflex**<sup>®</sup> bag technology.



**freeflex**<sup>®</sup>  
IV Bags & Solutions

All the features and advantages of **freeflex** – plus a unique Luer Lock injection port to facilitate needleless drug addition.



**Parenteral  
Nutrition**

Our parenteral nutrition products are used to help care for critically ill and chronically ill patients in hospitals, long-term care facilities, and at home.



To learn more about our expanding line of irrigation solutions, contact your IV Therapy Sales Representative or call Customer Service at **1.888.386.1300** | [www.FreseniusKabiIrrigation.com](http://www.FreseniusKabiIrrigation.com).